

## Chemical stability of the Varispenser® and the Top Buret

### Introduction

Bottle-top dispensers are used in the lab for dispensing a wide range of different solvents from glass or stainless steel containers. These dispensers have to meet various requirements. For example, they must not give off any substances which may disturb trace analysis, have cytotoxic properties, distort optical tests or influence chromatographic methods and residue analysis.



### Materials

Even after prolonged contact with the solvent, the materials of the dispenser must not be affected nor bind the solvent non-specifically. This means that there are very high demands on the chemical resistance of bottle-top dispensers.

The bottle-top dispensers Varispenser® and Varispenser® plus from Eppendorf are made of material which is particularly resistant to chemicals. Only parts made of PFA (perfluor alkoxy), PTFE (polytetrafluor ethylene), boron silicate glass 3.3

and platinum/iridium come into contact with the solvents. The adapter rings for the screw connection are made of PP (polypropylene).

The following pages include an example on the use of a Varispenser for residue analysis in the lab of a food manufacturer as well as a list of the materials of which the Varispenser is made and their chemical resistance.

### Dispensing ultrapure solvents for residue analysis at the Hipp plant in Pfaffenhofen, Germany

The residue analysis of foodstuffs places high demands on the solvents and the inertness of the lab equipment used to detect, for example, the smallest traces of pesticides.

Therefore, contamination caused by lab equipment has to be avoided under all circumstances.

A series of experiments were performed to test the suitability of the Varispenser for food analysis. Aim: to determine whether the Varispenser is inert to the solvents commonly used in this field.

### Two examples

#### Test A

100 ml of the solvent acetonitrile ( $\text{CH}_3\text{CN}$ ) were drawn from a larger supply with a Varispenser, evaporated in a rotation evaporator and then dissolved in 1 ml of the solvent i-octane. This sample was examined using gas Chromatography.

#### Result

No other peak could be found in

addition to the solvent peak.

Therefore, no substances had been released from the Varispenser.

#### Test B


A specific amount of the herbicide atrazine (50 pg/ml)<sup>1)</sup> was added to a sample treated as in test A as a comparative standard for pesticides. In order to detect minute amounts of pesticide, the background in the chromatogram must consistently be as low as possible.

#### Result

The atrazine peak can be clearly seen in the chromatogram without any ghost bands. This shows that minute amounts of pesticides can be detected with this system.

### Quotation from Mr. Lembacher, main lab residual analysis of the Hipp plant:

"The chromatograms recorded show that the dispenser is extremely inert. Fears about interference by plasticizers – even during long-term use – were totally unfounded."



<sup>1)</sup>(EG drinking water limit: 100 pg/ml)

#### Test A:

100 ml acetonitrile evaporated, dissolved in 1 ml i-octane, 0.5 µl injected.

#### Test B:

Like test A, except 50 pg/ml atrazine added as pesticide reference.

For each chemical, 2 figures are stated. The figure on the left is stability at a test temperature of +20 °C, the figure on the right is the stability at +50 °C. Salts were tested as almost saturated solutions. All data are recommendations without guarantee. 1 = resistant 2 = sensitive (raw material is affected after longer contact) 3 = incompatible

Materials	PP Adapter rings*				
	Top Buret/Varispenser can be used	Boron silicateglass 3.3	PTFE	PFA	
<b>Chemicals</b>					
<b>A</b>					
Acetaldehyde	1 1	1 1	1 1	1 1	1 3 <sup>*1</sup>
Acetic acid 50%	1 1	1 1	1 1	1 1	1 1
Acetone *5	1 1	1 1	1 1	1 1	1 1
Acetonitrile *5	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
Acrylonitrile	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
Adipic acid	1 1	1 1	1 1	1 1	1 1
Allyl alcohol	1 1	1 1	1 1	1 1	1 1
Aluminum chloride	1 1	1 1	1 1	1 1	1 1
Aluminum hydroxide	1 1	1 1	1 1	1 1	1 1
Amino acids	1 1	1 1	1 1	1 1	1 1
Ammonia	1 1	1 1	1 1	1 1	1 1
Ammonium chloride	1 1	1 1	1 1	1 1	1 1
Ammonium hydroxide 30%	1 1	1 1	1 1	1 1	1 1
n-Amyl acetate	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
Amyl alcohol	1 1	1 1	1 1	1 1	1 1
Amyl chloride	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Aniline	1 1	1 1	1 1	1 1	1 1
Aqua regia *2	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
<b>B</b>					
Barium chloride (BaCl <sub>2</sub> )	1 1	1 1	1 1	1 1	1 1
Benzaldehyde	1 2	1 1	1 1	1 2	1 1
Benzene	1 1	1 1	1 1	1 1	1 2 <sup>*1</sup>
Benzine	1 1	1 1	1 1	1 1	2 2 <sup>*1</sup>
Benzyl alcohol	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Biuret reagent	1 1	1 1	1 1	1 1	1 1
Boric acid	1 1	1 1	1 1	1 1	1 1
Bromine	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Bromoform (Tri Methanbromid)	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
n-Butanol	1 1	1 1	1 1	1 1	1 1
n-Butyl acetate	1 1	1 1	1 1	1 1	2 2 <sup>*1</sup>

Materials	PP Adapter rings*				
	Top Buret/Varispenser can be used	Boron silicateglass 3.3	PTFE	PFA	
<b>Chemicals</b>					
<b>C</b>					
Calcium chloride	1 1	1 1	1 1	1 1	1 1
Carbon disulphide	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Carbon tetrachloride	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Chloroacetic acid	1 1	1 1	1 1	1 1	1 1
Chromic acid 10%	1 1	1 1	1 1	1 1	1 1
Chromic acid 50% *2	1 1	1 1	1 1	1 1	2 2 <sup>*1</sup>
Chromic sulfuric acid, concentrated *2	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Cresol	1 1	1 1	1 1	1 1	1 2 <sup>*1</sup>
Cupric sulphate	1 1	1 1	1 1	1 1	1 1 <sup>*1</sup>
<b>D/E</b>					
Decahydronaphtalene	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Dibutyl phthalate	1 1	1 1	1 1	1 1	1 2 <sup>*1</sup>
Dichlorobenzene	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
Dichlorethane (Ethyl dichloride)*5	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
Dichlormethane (Methylene chloride) *5	1 1	1 1	1 1	1 1	3 3 <sup>*1</sup>
Diethylene glycol	1 1	1 1	1 1	1 1	1 1
Diethyl ether	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
Dimethylformamide	1 1	1 1	1 1	1 1	1 3 <sup>*1</sup>
1,4-Dioxan	1 1	1 1	1 1	1 1	2 2 <sup>*1</sup>
Ethanol 100% (Ethyl alcohol)	1 1	1 1	1 1	1 1	1 1
Ethyl acetate	1 1	1 1	1 1	1 1	1 1
Ethylene oxide	1 1	1 1	1 1	1 1	2 3 <sup>*1</sup>
<b>F/G</b>					
Formaldehyde 40%	1 1	1 1	1 1	1 1	1 1
Formic acid 98-100%	1 1	1 1	1 1	1 1	1 1
Fuel oil	1 1	1 1	1 1	1 1	1 1
Glycerol *4	1 1	1 1	1 1	1 1	1 1
Glycol	1 1	1 1	1 1	1 1	1 1

Materials	PP Adapter rings*				
	Top Buret/Varispenser can be used				Boron silicateglass 3.3
	PFA	PTFE			
<b>Chemicals</b>					
<b>H/I/J</b>					
n-Hexane	1 1	1 1	1 1	1 1	2 3 *1
Hydrochloric acid 35% *5	1 1	1 1	1 1	1 1	1 1
Hydrochloric acid 37% *5	1 1	1 1	1 1	1 1	1 3 *1
Hydrofluoric acid 40%	1 1	1 1	3 3	3 3	1 1
Hydrogen peroxide 30% *3	1 1	1 1	1 1	3 3	1 1
Iodine-potassium iodide sol.	1 1	1 1	1 1	1 1	1 1
Isobutanol (Isobutyl alcohol)	1 1	1 1	1 1	1 1	1 1
Isopropanol (Isopropyl alcohol)	1 1	1 1	1 1	1 1	1 1
Isopropyl benzene	1 1	1 1	1 1	1 1	2 3 *1
<b>L/M/N</b>					
Lactic acid/Lactate	1 1	1 1	1 1	1 1	1 1
Magnesium chloride (MgCl)	1 1	1 1	1 1	1 1	1 1
Mercury	1 1	1 1	1 1	1 1	1 1
Mercury (I) chloride	1 1	1 1	1 1	1 1	1 1
Methanol (Methyl alcohol) *5	1 1	1 1	1 1	1 1	1 1
Methyl propyl ketone	1 1	1 1	1 1	1 1	1 2 *1
Nitric acid 10%	1 1	1 1	1 1	1 1	1 1
Nitric acid 50%	1 1	1 1	1 1	1 1	2 3 *1
Nitric acid 70%	1 1	1 1	1 1	1 1	3 3 *1
Nitrobenzene	1 1	1 1	1 1	1 1	3 3 *1
<b>O/P</b>					
Octane/Iso octane	1 1	1 1	1 1	1 1	3 3 *1
Oil of turpentine	1 1	1 1	1 1	1 1	3 3 *1
Oxalic acid	1 1	1 1	1 1	1 1	1 1
Pentane (n-/Iso-) *5	1 1	1 1	1 1	1 1	3 3 *1
Perchloroethylene	1 1	1 1	1 1	1 1	3 3 *1
Perchloric acid 10%	1 2	1 1	1 1	1 1	1 3 *1
Phenol 100%	1 1	1 1	1 1	1 1	1 1
Phosphoric acid 85%	1 1	1 1	2 3	2 3	1 1
Potassium chloride	1 1	1 1	1 1	1 1	1 1
Potassium hydroxide 50%	1 1	1 1	1 2	1 2	1 1
Potassium permanganate	1 1	1 1	1 1	1 1	1 1
Propanol	1 1	1 1	1 1	1 1	1 1
Propylene glycol	1 1	1 1	1 1	1 1	1 1

\*1 PTFE adapter available

\*2 Pt-Ir can be easily loosened from the spring

\*3 Catalytic reaction with Pt-Ir spring

Materials	PP Adapter rings*				
	Top Buret/Varispenser can be used				Boron silicateglass 3.3
	PFA	PTFE			
<b>Chemicals</b>					
Propylene oxide	1 1	1 1	1 1	1 1	1 1
Pyridine	1 1	1 1	1 1	1 1	2 2 *1
<b>S</b>					
Salicylaldehyde	1 1	1 1	1 1	1 1	1 1
Salicylic acid	1 1	1 1	1 1	1 1	1 1
Scintillation cocktail	1 1	1 1	1 1	1 1	2 3 *1
Silver acetate	1 1	1 1	1 1	1 1	1 1
Silver nitrate	1 1	1 1	1 1	1 1	1 1
Sodium acetate	1 1	1 1	1 1	1 1	1 1
Sodium dichromate	1 1	1 1	1 1	1 1	1 1
Sodium hydroxide 50%	1 1	1 1	1 2	1 2	1 1
Sulphuric acid 60%	1 1	1 1	1 1	1 1	1 1
Sulphuric acid 98%	1 1	1 1	1 1	1 1	3 3 *1
<b>T/U</b>					
Tartaric acid	1 1	1 1	1 1	1 1	1 1
Tenside (Tween-, Triton X-, Brij-dilutions)	1 1	1 1	1 1	1 1	1 1
Tetrachloroethylene	1 1	1 1	1 1	1 1	3 3 *1
Tetrahydrofuran	1 1	1 1	1 1	1 1	1 1
Toluene	1 1	1 1	1 1	1 1	2 3 *1
Trichloroacetic acid 10%	1 1	1 1	1 1	1 1	1 1
Trichlorethane	1 1	1 1	1 1	1 1	3 3 *1
Trichlorethene	1 1	1 1	1 1	1 1	3 3 *1
Trichlorofluorethane *5	1 1	1 1	1 1	1 1	1 1
Trichloromethane (Chloroform)*5	1 1	1 1	1 1	1 2	2 2 *1
Triethylene glycol	1 1	1 1	1 1	1 1	1 1
Trifluoroacetic acid (fuming; strongest of halogenized acids)*5	3 3	1 1	3 3	3 3	3 3
Tripropylenglycol	1 1	1 1	1 1	1 1	1 1
Urea	1 1	1 1	1 1	1 1	1 1
<b>V/X</b>					
Vinylidene chloride	1 1	1 1	1 1	1 1	3 3 *1
Xylene	1 1	1 1	1 1	1 1	3 3 *1
<b>Z</b>					
Zinc chloride 10%	1 1	1 1	1 1	1 1	1 1
Zinc sulphate 10%	1 1	1 1	1 1	1 1	1 1

\*4 Liquid with high viscosity

\*5 Liquid with high vapor pressure;

gases leak (observe safety regulations)

## Material Varispenser/Top Buret

Part	Manual no.	Varispenser	Varispenser plus	Manual no.	Top Buret Bottletop Buret
		4960	4961		4965
<b>Direct contact to dispensing fluid</b>					
Valve block	15	PFA	PFA	9	PFA
Valve cock	11	-	PTFE	8	PTFE/PFA
Filling valve	17	ETFE	ETFE	10	ETFE
Discharge valve	11;12	ETFE	ETFE/PFA	-	ETFE/PFA
Discharge valve	10	PFA	PFA	7	PFA
Spring for valve	-	-	Pt-Ir	Pt-Ir	- Pt-Ir
Valve ball	16	Borosilicate (DURAN)	Borosilicate (DURAN)	-	Borosilicate (DURAN)
Cylinder	6	Borosilicate (DURAN)	Borosilicate (DURAN)	-	Borosilicate (DURAN)
Telescopic filling tube	-	FEP	FEP	12	FEP
<b>Indirect contact to dispensing fluid</b>					
Piston (2.5–10 ml)	5	PFA	PFA	-	PTFE/PFA
Piston (25–100 ml)	5	ETFE	ETFE	-	-
Piston holder	2	PP	PP	-	PTFE
Cylinder casing	3	PP	PP	-	PTFE/PFA
Protective cylinder sleeve	7	PTFE	PTFE	-	-
Valve block housing	14	PP	PP	2;10	PP
Discharge tube sleeve	9	PP	PP	5	PP
Discharge tube cap	13	PVDF	-	-	-
Air vent cup	18	PP	PP	11	PP
Volume adjustment knob	1	PP	PP	-	-
O-ring for valve cock protection	-	Viton	Viton	-	Viton
Volume setting knob	4	PP	PP	-	-
Discharge valve toggle	8	PP	PP	4	PP
Drying tube	-	PP	PP	-	PP
Wheel	-	-	-	3	PP
Display foil	-	-	-	1	PE

**DURAN** Borosilicate 3.3**ETFE** Tefzel ETFE (Ethylene tetrafluoroethylene)**FEP** Teflon FEP (Tetrafluoroethylene perfluoropropylene)**PFA** Teflon PFA (Perfluoro-alkoxy-PTFE-Copolymer)**PP** Polypropylene**PE** Polyethylene**PTFE** Polytetrafluoroethene**PVDF** Polyvinylidene fluoride**Pt-Ir** Platinum-Iridium

**eppendorf**  
In touch with life